

thermage® FLX

speed | efficacy | comfort



Product Presentation

A History of Innovation & Evolution



ThermaCool® TC system

- Pioneering non-invasive wrinkle reduction device

2002

Eye Tip 1.0 cm² (2002)
Face Tip 1.5 cm² (2004)
Face Tip 3.0 cm² (2005)
Eye Tip 0.25 cm² (2006)



Thermage® NXT system

- Second generation integrated system

2007

Body Tip 16.0 cm² (2008)



Thermage CPT® system

- Comfort pulse technology™

2009

Framed Tips (2009)
Total Tip 3.0 cm² (2013)



Thermage FLX™ system

- AccuREP™ technology

2017

Total Tip 4.0 cm² (2017)

The Latest Advancements



Faster
aLgorithm
eXperience



The Need For Speed



The Thermage FLX™ system has you covered. The **4.0 cm² Total Tip** available only on Thermage FLX can deliver a Thermage® treatment

25%
FASTER

compared to 3.0 cm² treatment tips. That is a **serious speed advantage!** What would you do with all that extra time?

Potential Value of Faster Treatments



**Faster treatment room
turnover**

**More revenue per
treatment hour**

More treatments per day


**Less time for provider
(per treatment)**

Less time for patient




OPTIMIZED ENERGY DELIVERY ALGORITHM

The Thermoage FLX™ system integrates tuning pulses into each individual radiofrequency energy pulse (REP). These tuning pulses automatically measure local impedance values and adjust energy output parameters to deliver consistent energy from one treatment site to the next.




Our Core Comfort Pulse Technology™

Intermittent RF pulses,
cooling bursts, and
vibration designed to
disrupt pain signals to
the brain



*Gate Control
Theory of Pain*



Intuitive Touchscreen Interface



One Handpiece For All Tips



- Lightweight
- Ergonomic
- Easy



We Can Treat the Eyes



Around the periorbital area and/or on the eye lid itself



We Can Treat the Face



Full face or
partial face
with a focus on
particular
problem areas

We Can Treat the Body

thermage[®]
FLX



With vibration
even when
using Body Tip
16.0

Environmentally Friendly Cryogen



Global Warming Potential (GWP)

1430  6



Example Scenarios



Thermage CPT/NXT

90 min/Tx
* 2 Tx/wk
* 4 wk/mth
= 12 Tx hrs/mth

8 Tx/mth
* Avg \$2,500/Tx
= \$20,000/mth

Thermage FLX

45 min/Tx
* 2 Tx/wk
* 4 wk/mth
= 6 Tx hrs/mth

8 Tx/mth
* Avg \$2,500/Tx
= \$20,000/mth
+ extra 6 hrs Tx
room time/mth

Avg Tx room revenue = \$1,000/hr

6 Tx hrs/mth * \$1,000/hr * 12 mths =

\$72,000 per year

Example scenarios provide estimates and general guidance to help you determine whether you desire to further investigate the benefits of purchasing and using selected Solta products. EXAMPLE SCENARIOS ARE NOT A GUARANTEE OF ACTUAL FINANCIAL PERFORMANCE.

INDICATION

The **radiofrequency energy only delivery components of the Thermage FLX™ system** and accessories are indicated for use in:

- Non-invasive treatment of periorbital wrinkles and rhytids including upper and lower eyelids; and
- Non-invasive treatment of wrinkles and rhytids.

The **simultaneous application of radiofrequency energy and skin vibration by the Thermage FLX system** and accessories are indicated for use in:

- Non-invasive treatment of periorbital wrinkles and rhytids;
- Non-invasive treatment of wrinkles and rhytids; and
- Temporary improvement in the appearance of cellulite.

IMPORTANT SAFETY INFORMATION

- Do not use the Thermage FLX system in patients with cardiac pacemakers or other active implants.
- Improper use of the Thermage FLX system may cause personal injury or damage to the system. See the Technical User's Manual for detailed directions, proper use, and full risk and safety information.